# MODE D'EMPLOI - INSTRUCTIONS FOR USE

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*Translation into English from original in French*

## D. INSTRUCTIONS FOR USE

### 1. SAFETY INSTRUCTIONS

#### 1.1. USE

Duplicator for gel DUCAGEL Z1 has to be operated only by dental professionals.

DUCAGEL Z1 has to use in the conditions and for the uses described in this manual. Any other use that will be made of it (such as firing other materials, heating food products, etc.) is considered inappropriate. In this case, the manufacturer declines all liability and immediately cancels the equipment warranty.

According to the product’s instructions and functionalities, **it is imperative to read carefully the instructions for use before installation and use.** Retain the manual for future reference.

#### 1.2. SAFETY INSTRUCTIONS

Pictograms’ meaning used in this instructions for use:

 General warning	 Electrical hazard	 Burn hazard - Hot surface
---	---	---

Observe following safety instructions:

	<ul style="list-style-type: none"> <li>- During operation and after use, <b>the surface may reach high temperatures: do not touch</b></li> <li>- Liquefied gel may reach high temperatures (in function to settings by the user): do not touch when exit by the spout</li> </ul>
	<ul style="list-style-type: none"> <li>- The equipment has to be connected to an <b>electrical installation in accordance</b> with the applicable standard in the country in which it is being used. This installation has to provide <b>protection against overcurrent, overload and ground faults</b></li> <li>- <b>Connect it directly to the mains socket.</b> Do not use multiple sockets</li> <li>- Ensure that the <b>connection cables</b> are in perfect condition to prevent short circuits</li> <li>- Before any <b>technical intervention</b>, the equipment has to be switched off (with the breaker) and disconnected from the mains power supply</li> <li>- <b>Do not spill liquid</b> on the ventilation openings or inside the equipment</li> </ul>
	<ul style="list-style-type: none"> <li>- The equipment must be placed on a <b>flat surface</b>, keeping sufficient clearance around the unit to <b>ensure proper ventilation</b></li> <li>- <b>Do not block the ventilation openings</b> to avoid overheating</li> </ul>
	<ul style="list-style-type: none"> <li>- <u>Insulation plate</u> placed under tank contains <b>insulating material</b> composed of high temperature insulation wool (alkaline earth silicate), a substance not classified as dangerous according to Regulation 1272/2008/EC (CLP). Handle carefully</li> <li>- If handling, dust can be released: remove with a vacuum cleaner - do not blow or use compressed air</li> </ul>

- Use only **genuine spare parts**. The use of non-original spare parts voids any warranty for your device
- **Do not introduce stranger objects inside** the equipment during the maintenance operations. It is forbidden to modify the material without authorization
- **Do not execute different operations of maintenance** from those brought back in the manual. Whichever operation not included in this manual, can involve risks
- For whichever information concerning installation, maintenance and use, contact the **customer's service UGIN DENTAIRE**

### 1.3. WASTE DISPOSAL

Observe the waste disposal regulations for electrical and electronic equipment: do not dispose of with household waste.

The unit is subject to Directive 2012/19/EU on waste electrical and electronic equipment and to the laws of the country in which it is being used, and must be disposed of in accordance with applicable regulations.

Insulation plate is composed of high temperature insulation wool (alkaline earth silicate), a substance not classified as dangerous according to Regulation 1272/2008/EC (CLP), and classified as non-hazardous waste for disposal. However, such a waste is normally dusty and so must be properly packaged before for disposal. Check any national or regional applicable regulations.

## 2. EC DECLARATION OF CONFORMITY

The equipment's EC declaration of conformity is provided as appendix of instructions for use.

### Content of the declaration EC:

- Name and address of the manufacturer:  
UGIN DENTAIRE. 25 rue de la Tuilerie. 38170 SEYSSINET-PARISSET – FRANCE
- UGIN DENTAIRE declares, on its own exclusive responsibility, that the product:
  - Duplicator for gel DUCAGEL Z1
  - Serial number: *see the equipment's EC declaration of conformity*
- Fulfills the essential requirements of health and safety of the following directives:
  - Machinery directive 2006/42/EC
  - Electromagnetic compatibility directive 2014/30/EU

### 3. TECHNICAL INFORMATION

#### 3.1. TECHNICAL DATA

Duplicator for gel DUCAGEL Z1 has to be operated only by dental professionals and in the conditions described in this manual (see chapter 1.1 and 3.2).

<b>TECHNICAL DATA</b>	
<b>Maximum vat temperature</b>	99°C
<b>Mixing speed</b>	30 revolutions/minute and includes electronic torque control to prevent damage to the drive unit
<b>10 programs</b>	- <b>No 0</b> : continuous mixing action during the low-temperature holding level - <b>No. 1 to No. 9</b> : intermittent mixing action
<b>Delayed start time up</b>	to 99 hours
<b>Voltage</b>	230 V – 50/60 Hz
<b>Power consumption</b>	850 W
<b>Dimensions</b>	
- height (lid closed)	510 mm
- height (lid open)	820 mm
- width	325 mm
- depth	425 mm
<b>Vat dimensions</b>	
- internal diameter	240 mm
- height	220 mm
- capacity	6 kg of gel
<b>Weight</b>	23.5 kg
<b>Noise level</b>	< 70 dB(A)

### 3.2. CONDITIONS OF USE, TRANSPORT AND STORAGE

#### ▪ **Authorised conditions of use**

- Ambient temperature: +5°C to +40°C
- Humidity range: maximum relative humidity 80% for temperatures until 31°C, without condensation and decrease linear until 50% to 40°C, without condensation
- Ambient pressure: up to an altitude of 2000 m above sea level

#### ▪ **Authorised conditions of storage**

- Ambient temperature: -20°C to +65°C
- Humidity range: maximum relative humidity 80%
- Ambient pressure: up to an altitude of 2000 m above sea level

#### ▪ **Conditions to transport**

- Wait the complete cooling before packing/transporting
- Use the original packaging only, with the protective elements
- Take hold the equipment by its base
- Install the equipment in the packaging with the protective elements.

## 4. INSTALLATION AND START-UP

### 4.1. UNPACK

- Unpack DUCAGEL Z1 and check that it is in perfect order
- Any defects can be reported to the carrier
- Remember to get the carrier to sign the delivery note

*NOTE: keep the original packaging to transport the equipment (see chapter 3.2)*

### 4.2. INSTALLATION

- Install in a room aired. Observe the conditions indicated in chapter 3.2
- Place the equipment on a level and aired surface. Keep a distance from other objects around not less than 10 cm. Keep it away from heat sources (for example, radiators and/or other equipment that release heat). Avoid vibrations and shocks
- Do not install in a place with explosion hazard
- Its forbidden placing inflammables, toxic, volatile or explosives around the equipment

**Grounding is mandatory** to avoid accidents and the dangers of overheating the power line.

Due to the power of equipment: **connect it directly to the mains**, do not use multiple sockets (risk of electrical interference)

### 4.3. START-UP

NOTE: for marks, see chapter A – equipment description

- Connect the power cord to the power socket (at the back panel)
- Make sure that the mains voltage complies with the voltage indicated on the rating plate, connect the power cord to a 10/16 A grounded power socket
- Fill the tank with gel (max. 6 kg)
- Activate the ON/OFF switch (at the back panel)
- The temperature in the tank and the program number will appear after a few seconds
- The programs can now be configured (see chapter 5)

## 5. CONFIGURATION AND USE

NOTE: for marks, see chapter A – equipment description

### 5.1. CONFIGURATION

For select the program: use the keys  and 

#### 5.1.1. DELAYED CYCLE START

- Press the key  to start the delayed cycle start configuration
- The letter **d** appear on the program number display (n° 1)
- The number **0.0** appear on the temperature/time display (n° 2)
- The time adjustment range extends up to 00 hours 00 seconds to **99 hours and 00 seconds**
- To be efficient, **the value has to more 60 seconds**

NOTE: number displayed (0.0) is made up of two parts

0.0
hours   seconds

- Select the time value with the keys  and 
- As soon as pressing the key  :
  - if the programmed settings is 0.0 : the cycle starting
  - for a delayed cycle start: the countdown starting (before the cycle required will start)

## 5.1.2. PARAMETERS SETTING

### 5.1.2.1. LOW STABILISATION TEMPERATURE

- Press the key  to set the low stabilisation temperature: the letter **b** appear to the left on temperature/time display (No. 2)
- The adjustment range extends up to **30°C to 75°C (preset to 45°C)**
- Select the value with the keys  and 

### 5.1.2.2. HIGH STABILISATION TEMPERATURE

- Press again the key  to set the high stabilisation temperature: the letter **h** appear to the left on temperature/time display (No. 2)
- The adjustment range extends up to **75°C to 99°C (preset to 94°C)**
- Select the value with the keys  and 

### 5.1.2.3. STABILISATION TIME OF THE TOP LEVEL

- Press again the key  to set the stabilisation time of the top level: the letter **t** appear to the left on temperature/time display (No. 2)
- The adjustment range extends up to **10 minutes to 30 minutes maximum (preset to 10 minutes)**
- Select the value with the keys  and 

### 5.1.2.4. VALIDATION

- Press the key  to validate all parameters set: low and high stabilisation temperature and stabilisation time of the top level
- Repeat the operation for **each programs**

## 5.2. USE

- Select the program with the keys  and 
- Press the key  to **start the cycle**
- **Start to heating process:** the heater indicator lights up (*No. 11*)
- As soon as the temperature in the tank reaches 35°C, the **intermittent rotation of the knife** starts and the drive unit indicator lights up (*n° 9*)
- The **temperature remains constant** as soon as the high temperature set point is reached, until the end of the stabilisation time of the top level (t)
- At the end of the stabilization time of the top level, the **fan starts** and its indicator lights up (*No. 10*): the mixing continues but the heating is inactive
- Once the low temperature set point (b) has been reached, the fan stops and the motor runs only intermittently during the low stabilization time (except for program No. 0)

**CAUTION: for safety reasons, THE CYCLE NOT START IF THE LID IS OPEN**

In this case, *OPEN* appear on the display: the motor and the heater process are inactive

- **NOTE**
  - **Stabilisation phase:** indicated by a dot at the bottom right of the program number
  - **Program No. 0:** the knife rotate continuously until the end of the low holding level
  - If necessary, you can **switch directly to the low-temperature stabilisation phase (b):** press the key  then the key 
  - To **return to the beginning of the cycle** press the key  then 
  - To press the key  **interrupts the current program:** stoppage drive unit, fan and heating cycle, and the temperature drop naturally
- **RECOMMENDATION**
  - **Break the gel down into small pieces** so that the drive unit can restart more efficiently and to reduce the time required to reach liquefaction temperature
  - **After a long period of stabilisation,** a small amount of water should be added to compensate for any natural evaporation

## 6. MAINTENANCE AND TROUBLE-SHOOTING

**Before to carry out the maintenance work and/or cleaning:**

- **disconnect from the main power supply** (*switch OFF and unplug to mains socket*)
- **wait the cooling at room temperature** after use (*risk of burns and/or damage*)

### 6.1. MAINTENANCE

Observe the following recommendations to ensure the longevity of the device:

- Clean the machine outside with a soft, dry cloth
- Do not use abrasive, solvents or flammables products for cleaning (outside and inside)
- Periodically, **clean the inside of the tank with water to remove any residue gel**. If necessary, remove the knife (**warning: very sharp blades**)

### 6.2. TROUBLE-SHOOTING

<i>Operational failure</i>	<i>Cause</i>	<i>Procedure to be followed</i>
<b>Equipment does not heat</b>	<b>1. Program setting (temperature)</b>	- Check the program's parameters setting: low and high stabilisation temperature value (see chapter 5.1.2)
	<b>2. Resistance / electronic card</b>	- Resistance or electronic card is defective or improperly connected : contact the after-sales service

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