**MODE D'EMPLOI - INSTRUCTIONS FOR USE**

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*Translation into English from original in French*

## D. INSTRUCTIONS FOR USE

### 1. SAFETY INSTRUCTIONS

#### 1.1. USE

**INTEGRAL 2 is intended for carry out an electrolytic polishing by electro-substraction.**

**This equipment has to be operated only by dental professionals.**

**INTEGRAL 2 has to use in the conditions and for the uses described in this manual.** Any other use that will be made of it is considered inappropriate. In this case, the manufacturer declines all liability and immediately cancels the equipment warranty.

**It is imperative to carefully read the instructions for use before installation and use.** Retain the manual for future reference.

#### 1.2. SAFETY INSTRUCTIONS

Pictograms' meaning used in these instructions for use:

	General warning		Electrical hazard
---	-----------------	---	-------------------

**Observe following safety instructions:**

	<ul style="list-style-type: none"> <li>- The unit has to be <b>connected to an electrical installation in accordance</b> with the applicable standard in the country in which it is being used. This installation has to provide protection against overcurrent, overload and ground faults.</li> <li>- <b>Connect it directly to socket.</b> Do not use multiple sockets.</li> <li>- Ensure that the <b>electric cord</b> is in perfect condition to prevent short circuits.</li> <li>- Before any <b>technical intervention</b>, the unit has to be switched off and disconnected from the mains power supply.</li> <li>- <b>Do not spill liquid</b> inside the equipment.</li> </ul>
	<ul style="list-style-type: none"> <li>- Not exceed the maximum tank's volume (2 litres).</li> </ul>
	<ul style="list-style-type: none"> <li>- Use only <b>genuine spare parts</b>. The use of non-original spare parts voids any warranty for your device.</li> <li>- <b>Do not introduce stranger objects inside</b> the unit during the maintenance operations. It is forbidden to modify the material without authorization.</li> <li>- <b>Do not execute different operations of maintenance</b> from those brought back in the manual. Whichever operation not included in this manual, can involve risks.</li> <li>- For whichever information concerning installation, maintenance and use, contact the <b>customer's service UGIN DENTAIRE</b>.</li> </ul>

#### 1.3. WASTE DISPOSAL

Observe the waste disposal regulations for electrical and electronic equipment: do not dispose of with household waste.

The unit is subject to Directive 2012/19/EU on waste electrical and electronic equipment and to the laws of the country in which it is being used, and must be disposed of in accordance with applicable regulations.

## 2. EC DECLARATION OF CONFORMITY

The equipment's EC declaration of conformity is provided as appendix of instructions for use.

### Content of the declaration EC:

- Name and address of the manufacturer:  
UGIN DENTAIRE  
25 rue de la Tuilerie  
38170 SEYSSINET-PARISSET – FRANCE
- UGIN DENTAIRE declares, on its own exclusive responsibility, that the product:
  - Electrolytic polishing **INTEGRAL 2**
  - Serial number: *see the equipment's EC declaration of conformity*
- Fulfills the essential requirements of health and safety of the following directives:
  - Machinery directive 2006/42/EC
  - Electromagnetic compatibility directive 2014/30/EU

## 3. TECHNICAL INFORMATION

**INTEGRAL 2** is intended for carry out an electrolytic polishing to by electro-substraction. This equipment has to be operated only by dental professionals and in the conditions described in this manual.

### ▪ ELECTROLYTIC POLISHING

The **INTEGRAL 2** generator uses the principle of electro-substraction: the workpiece is submerged in an acid bath, the high current flowing through it allows a quick and effective polishing (5 to 15 minutes).

To achieve a polished surface, electric abrasion decreases haight a in a process known as electrolytic polishing or electropolishing. This involves electrolysis of the workpiece, attacking the metal, which is then covered by a thin liquid coating with a high metallic ion content. As these ions are less mobile than hydrogen ions, this coating enhances the resistance of the electrolyte to the electric current passing through it. The surface coating varies in thickness, thickest in the hollows and thinnest at the peaks. This makes for a higher current density at the peaks, which thus dissolve more quickly.

<b>TECHNICAL DATA</b>	
<b>Dimensions</b>	
- Width	330 mm
- Height	175 mm
- Depth	215 mm
<b>Weight (empty)</b>	6 Kg
<b>Supply voltage</b>	230 V
<b>Frequency</b>	50 / 60 Hz
<b>Secondary voltage</b>	18 V
<b>Maximum amperage</b>	10 A
<b>Timer</b>	0 to 30 minutes
<b>Circuit-breaker</b>	10 A at secondary
<b>Fuse</b>	2 A at primary
<b>Composition tank</b>	PVC
<b>Capacity tank</b>	2 liters
<b>Noise level</b>	< 70 db(A)

### 3.1. ACCESSORIES

- **Accessories delivered** (included in the packaging):

<b>DESCRIPTION</b>	<b>CODE UGIN</b>
Copper wire (1 meter)	FCE05DIV0001

- **Accessories recommended** (not delivered – to order separately):

<b>DESCRIPTION</b>	<b>CODE UGIN</b>
- Electrolytic liquid – 1L	PLS03ALR0001
- Electrolytic liquid – 3L	PLS03ALR0002

## 4. INSTALLATION AND START-UP

### 4.1. UNPACK

- Unpack **INTEGRAL 2** and check that it is in perfect order.
- Any defects can be reported to the carrier.
- Remember to get the carrier to sign the delivery note.

*NOTE: keep the original packaging to transport.*

## 4.2. INSTALLATION AND USE

- a) Install in a room aired, and place INTEGRAL 2 on a level surface.
- b) Fill the tank with the polishing liquid.

**IMPORTANT: never exceed the maximum level (max 2L)**

We recommend using the **UGIN DENTAIRE polishing bath: electrolytic liquid** (see chapter 3.1). It is exceptionally effective because it heats up rapidly. In addition, it has a particularly good stability over time.

**IMPORTANT: as the polishing quality depends primarily on the kind of solution used. We disclaim results obtained with the solutions other than our own.**

- c) Connect the workpieces to polishing  
A copper wire, supplied with the device, is used to connect the copper bar to the workpiece. Check that proper contact has been made between the two.
  - The copper bar should be cleaned regularly with fine-grained emery paper.
  - When the copper wire becomes too thin, it should be replaced.
- d) Switch on.
  - Make sure that the mains voltage complies with the voltage indicated on the rating plate, connect the power cord to a 10/16 A grounded power socket.
  - Press the power switch ON/OFF (rear side).
- e) Turns the timer (0-30 minutes) to the required time.  
To achieve a good polish, the average operating time is 10 to 15 min.

**CAUTION: do not manually reset the timer to zero** (risk of damage to the equipment).  
If necessary, switch off with the power switch ON/OFF.

- f) Adjust the intensity with the power control button.  
The current densities most frequently used in electropolishing are around 0,5 A per cm<sup>2</sup>.  
The surface area of frameworks may vary considerably, but on an average they are 6 to 10 cm<sup>2</sup>. Hence the amperage for a given framework is to 3 to 5 A.  
For our machine, we recommend **limiting the total area to be polished to 20 cm<sup>2</sup>, or two to three frameworks, for best results.**
- g) When the polishing time expired, recover the workpieces in the tank.  
Use appropriate gloves (observe the recommendations of the manufacturer of polishing liquid - see technical data sheet and/or safety data sheet)

### **SPECIAL CARE MUST BE TAKEN WITH THE CLASPS**

Excessive current density or overly long exposure time may dangerously narrow the cross-section of such protruding parts. At the back of the palate, which is by definition concave, it is normal to observe a lower polish.

**A little mechanical polishing may be required to achieve the desired result.**

## 5. MAINTENANCE

**Before to carry out the maintenance work and/or cleaning: disconnect from the main power supply (switch OFF and unplug to mains socket).**

Observe the following recommendations to ensure the longevity of the device:

- Do not use abrasive, solvents or flammables products for cleaning (outside and inside).
- Clean the machine outside with a soft, dry cloth.
- **Change the polishing liquid periodically.**  
*NOTE: For periodicity, observe the recommendations of the manufacturer of liquid.*
- The **copper bar should be cleaned regularly** with fine-grained emery paper.
- When the **copper wire** becomes too thin, it should be **replaced**.

### 5.1. TROUBLE-SHOOTING

<b>PROBLEM</b>	<b>PROCEDURE TO FOLLOWED</b>
<p><b>SHORT-CIRCUIT</b> If the workpiece comes in contact with the annular electrode for several seconds, the thermal circuit-breaker breaks the current.</p>	<p>Should this occur, wait for one minute and then reset the circuit-breaker. The device is now ready for operating again.</p>
<p><b>POOR CONDUCTION</b> If the ammeter shows no change, despite the rotation of the potentiometer to the right.</p>	<ul style="list-style-type: none"> <li>- Check the connection between the bar and the workpiece.</li> <li>- Check the 10 A circuit-breaker.</li> <li>- Check the fuse.</li> </ul>
<p><b>POLARITY REVERSAL</b> If the workpiece builds up a copper layer:</p>	<p>Correct the polarity to eliminate the coating.</p>

*If the problem persists, contact the after-sale service*

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