# DUCATRON QUATTRO

CE

**UGIN'**  
dentaire

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# USER GUIDE

<b>INTRODUCTION</b>	<b>13</b>
<b>TECHNICAL CHARACTERISTICS</b>	<b>14</b>
<b>INSTALLATION AND START-UP</b>	<b>15</b>
<b>CONTROL PANEL DESCRIPTION</b>	<b>16</b>
<b>OPERATIONAL TEST WHEN EMPTY</b>	<b>17</b>
<b>ARM BALANCING</b>	<b>17</b>
<b>ADVICE ON USING CRUCIBLES</b>	<b>18</b>
<b>CYCLE SEQUENCE</b>	<b>19</b>
<b>OPERATIONAL FAULTS</b>	<b>20</b>
<b>SERVICING AND PRECAUTIONS FOR USE</b>	<b>21</b>

*Translation into English from original in French*

## INTRODUCTION

This new generation induction centrifuge has been designed better to respond to your needs.

The direct movement transmission ensures constant start power along with a very low sound level during injection.

The productivity reached by its new induction generator is higher than that previously generated on the Series 3 DUCATRON.

Particular attention has been given to safety, accessibility and design.

Several sensors check that the machine is operating correctly, as follows:

- A warning light on the control panel indicates any lack of water in the cooling circuit, to protect the self-induction system.
- Then the light is on the heating system becomes inoperable.
- A warning light on the front face of the machine indicates any generator malfunction.

Two other LEDs are part of the control system:

- One confirms the arm position (thus authorizing coil movement)
- and the other confirms the cap as being in position, thus allowing injection to take place in complete safety.

Another feature is that the DUCATRON Quattro compact dimensions ensure easy installation in all types of dental prosthesis laboratories.

The instrument design gives easy access to the various control elements.

Only the parts at the front of the machine are able to be removed, thus rendering any intervention simple and rapid.

# TECHNICAL CHARACTERISTICS

## OVERALL DIMENSIONS:

- WIDTH:	600 mm
- HEIGHT:	980 mm
- DEPTH:	620 mm
- WEIGHT:	95 Kgs

## ELECTRIC POWER SUPPLY:

- TOTAL POWER:	2500 W
- ABSORBED GENERATOR POWER:	2000 W
- ELECTRIC POWER VOLTAGE:	230 V
- FREQUENCY:	50 Hz

## CHARACTERISTICS:

- GENERATOR FREQUENCY: 130 KHz.
- FUSION CAPACITY: 60g precious or non-precious alloy.
- MAX TEMPERATURE: 1500/1600°c
- ADJUSTABLE POWER INPUT: from 20 to 100%
- CENTRIFUGE SPEED: 450 rpm
- TWO TYPES OF STARTING MODE:
  - . "Moderate" for precious alloys (cylinders 1X and 3X).
  - . "Rapid" for non-precious, semi-precious and precious alloys (cylinders 6X and 9X).
- MEAN FUSION TIME FOR 25g NEW ALLOY: 1 minute
- WATER LEVEL WARNING LIGHT (Cooling circuit).
- GENERATOR MALFUNCTION WARNING LIGHT.
- ON AUTOMATIC STANDBY: After 10 minutes without control buttons being touched
- CAPACITIVE TOUCH CONTROL PANEL (Light touch).
- CLOSED CIRCUIT INDUCTOR COOLING SYSTEM:
  - water pump, filling 20 litre minimum to 22 litre and water sensor.
- AUTOMATIC ARM RETURN TO INITIAL POSITION AT CYCLE END.

## ACCESSORIES:

- 5 CRUCIBLES: 2 GREY (Silicon carbide) and 3 WHITE (Silica).
- 1 VITRIFIED CARBON SLEEVE.
- 1 QUARTZ BATON.
- 6 CRADLES (1X-3X-6X-9X-D42-D92).
- 4 METAL CYLINDER ASSEMBLIES (1X-3X-6X-9X).
- 1 GRIPPER FORK.
- 1 SET OF BREATHER CONNECTORS.
- 4 ADJUSTABLE FEET (For machine stability).


# INSTALLATION AND START-UP

## **PLEASE NOTE!**

The DUCATRON Quattro is designed to be used only in dental prosthesis laboratories. Any other type of use that the machine might be subjected to immediately invalidates the manufacturer's guarantee.

1. When unpacking, check machine condition so as to inform the transporter and supplier of any damage or deterioration to the machine or its accessories. Do not forget to sign the transport form on which any reservations can be made, and confirm these to the transporter within 48 hours by recorded delivery letter with acknowledgement of receipt.
2. After fitting the four adjustable feet in place on the underside of the chassis, position the DUCATRON QUATTRO on a flat surface.  
Allow some 5 to 10 cm space behind the machine.  
Use foot adjusting screws to level and stabilize machine.  
Please note: This operation is important: any instability brings about vibrations and considerably increases the machine sound level during any centrifuge operation.
3. Take out plugs, and then unscrew the four machine panel screws using a crosshead screwdriver.  
Remove panel and disconnect earth lead (yellow and green).
4. Swing out the reservoir holding bracket by partially loosening the butterfly nut.
5. Free the band holding the cover in place so as to gain access to container.
6. Fill reservoir with 20 liters of water minimum to 22L maximum.
7. Put container - then pump and cover - back in place.
8. Close up holding band while ensuring all elements are correctly positioned.
9. Connect up panel earth lead.
10. Close cover, and then insert the four plastic plugs.
11. Open upper cover and take out arm packing foam.
12. Having checked that the mains voltage is that specified on the machine identification label, connect power cable to a standard 10/16 A power plug that has an earth prong.
13. Throw bipolar isolating switch - to the right on machine upright - to ON.
14. Push ON button on control panel.  
The cooling pump starts.
15. The DUCATRON Quattro is now ready to operate.

# CONTROL PANEL DESCRIPTION

- 1 – "ON" Power button
- 2 – Power on warning light.
- 3 – Water level warning light (Lights up when water low).
- 4 – Generator warning light (Lights up in case of generator plate failure).
- 5 – "Moderate" start warning light.
- 6 – Start mode selector switch.
- 7 – "Rapid" start warning light.
- 8 – Generator power level warning lights (from 20 to 100%).
- 9 – Power selection switch.
- 10 – "FUSION" warning light
- 11 – "INJECTION" warning light
- 12 – "STOP" button (cycle stop or phase under way stop).
- 13 – "INJECT" button (coil descent control and arm rotation start).
- 14 – "FUSION" (alloy fusion start control).
- 15 – "SELF" [inductor coil] (inductor engage control when arm is in median position).
- 16 – Arm position indicator (Lit in median position).
- 17 – Cover locked warning light.

## OPERATIONAL TEST WHEN EMPTY

To familiarize you with the DUCATRON Quattro, we advise you to carry out your first test "on empty" by following the directions below.

1. Place the 30mm Ø stainless steel cylinder in the smallest of the cradles delivered as part of the accessories.
2. Place the cradle + cylinder assembly on the machine arm.
3. Position the silica crucible fitted with its vitrified carbon sleeve in its designed place on the arm.
4. Switch on the power to the machine using the "**ON**" switch on the control panel.
5. Bring the arm to its median. A signal sounds and a warning light confirms correct positioning.
6. Activate "**SELF**" switch.
7. Close cover.
8. After setting generator power to 80% (4 lights illuminated), operate "**FUSION**" switch.
9. Wait about 15 seconds (carbon sleeve goes red) then operate "**INJECT**" switch.
10. After some fifteen seconds push on "**STOP**" button to stop cycle.
11. The cover automatically unlocks some seconds later.

## ARM BALANCING

The DUCATRON Quattro arm balancing system is identical to that well-tried on the DUCATRON Series 3 machine.

It is effective because of the six cradle models designed for the various cylinder types.

As a function of cradle size, a specific counterweight is fitted to maintain a total constant weight (1150g ±70g).

For information the table below shows the theoretical weights of the different (clean) cylinders along with the weights of the corresponding cradles.

CYLINDER DIMENSIONS IN MM		CYLINDER WEIGHT IN GRAMS		CRADLE WEIGHT IN GRAMS
Cylinder dia	Height	Metal	Non-metal	
30	55	100	54	1110
40	55	170	96	1035
50	55	230	150	979
65	55	355	255	835
80	55	500	385	700
90	70		623	539

### PLEASE NOTE!

Never start to rotate arm without having checked for the presence of a cradle and its cylinder. Having no cradle and cylinder on the arm will create strong vibrations in the machine and could damage internal parts.

# USING CRUCIBLES

The DUCATRON Quattro is delivered with two types of crucible and a vitrified carbon sleeve.

– SILICON CARBIDE CRUCIBLE (grey)

It can be used for non-precious metals.

Its composition means it can be used up to temperatures of 1500°C to 1600°C.

The silicon carbide helps prevent the tongue formed during injection between crucible bottom and lip sticking to the crucible.

– SILICA CRUCIBLE (white)

It can be used with precious and non-precious alloys. Its maximum operating temperature (1400°C à 1500°C) is less than that for silicon carbide crucibles.

On the other hand, its cost price is lower.

– VITRIFIED CARBON SLEEVE (or vitreous carbon)

Designed for using precious or semi-precious alloys, it can be used with of silica crucible.

Being smooth it protects the alloy from any contamination by cup components.


**PLEASE NOTE!**

Never use the vitrified carbon sleeve for fusing non-precious alloys.

**ADVICE ON USING CRUCIBLES:**

Casting incidents are generally linked to the choice or condition of crucible used, in terms of possible fissures, distortion or wear. To avoid such problems, it is recommended to respect some simple rules:


1. Check on the type of crucible recommended by the alloy manufacturer. If no particular is detailed, use cups mentioned above.
2. If the crucible has already been used, check its condition.
3. Watch over the fusion process constantly to avoid any alloy overheating.  
Any prolonged holding of the fusion temperature brings about distortion and premature cup wear (Fig. 1).

4. Never allow cast metal that has not been injected to cool down in the crucible.
5. Do not use feeders that are too big, and at filling avoid bits of alloy jamming between the crucible walls (Fig. 2).

Please note: It is also recommended that you first put feeder pieces, then new alloy slugs, into the crucible.

### FEEDER PIECES

## **CYCLE SEQUENCE**

- Flip isolating switch to "ON".
- Push ON button on control panel.
- Use button "**C**" to select type of arm start mode.
- Load crucible with the desired amount of alloy (60g maximum).
- Insert crucible into its housing on the arm.
- Check arm centering. A signal sounds to confirm correct positioning (Light on front panel also lights up).

**Please note:** When the sound signal and the warning light are not activated, the "**SELF**" [coil] button remains inoperative

- Use the gripper fork to place the cylinder and cradle assembly in the arm compartments.
- Activate the fusion coil movement ("**SELF**" button). Two "beeps" confirm the inductor position.
- Parameterize generator power ("**P**" button from 1 to 5 warning lights lit for max power)
- Close hinged cover.
- Start fusion ("**FUSION**" button).
- Slugs melt into each other then form a homogenous ball.
- Injection temperature is then reached.
- Press "**INJECT**", cover locks (warning light lit), inductor lowers to release crucible and centrifuge starts.
- Press "**STOP**" to stop rotation after 10 to 15 seconds.
- Wait for the cover safety lock to release cover (about 20 seconds).
- Open cover and use gripper fork to recover cylinder and cradle assembly.
- Use anti-heat gloves or a suitable gripper to take out crucible.

## OPERATIONAL FAULTS

### WATER FAULT:

Water fault warning light (  ) is on.

Water volume to cool inductor is insufficient.  
Fusion is de-activated in this case to protect coil.

Filling reservoir: See section on servicing and precautions on use

If fault persists contact the after-sales service UGIN'DENTAIRE.

### GENERATOR FAULT:

Generator fault warning light (  ) is on.

Contact the after-sales service UGIN'DENTAIRE.

# SERVICING AND PRECAUTIONS FOR USE

## SERVICING:

Except for cleaning the machine outside using products without solvents, the DUCATRON Quattro needs no particular servicing operation.

However it is recommended that the cooling circuit reservoir water level is checked about every three months.

### To carry out this operation:

- Cut off power to the machine (isolating switch on right to "OFF").
- Take out front lower panel without forgetting to disconnect earth lead.
- Undo band that secures cover assembly.
- Release cover to gain access to reservoir opening.
- Add any needed water so that the water level is about 5 cm below the upper edge.
- Put back the various parts, checking that the cover band is correctly positioned and locked up.
- Re-connect panel earth lead.
- Close up machine (Four quarter-turn screws + four plugs).
- Reset mains power.
- The DUCATRON QUATTRO is operational again.

## PRECAUTIONS FOR USE:

Before any use, we would advise you to check that no object or residue is present in the machine arm.

It is essential that the cradle fits perfectly with the cast cylinder for any operation to be good. It is best to carry out this check when the machine is cold.

The machine is delivered with four cylinder holder cradles (DEGUSSA Standard) and two others that can be used with 40mm Ø x 55mm long and 90mm Ø x 70mm long cylinders.

Other types can be used on condition that their dimensions and weights remain compatible with the accessories of the machine.

To prevent any risk of alloy projection during centrifugation, several points should be checked:

- Distance between cup lip and cylinder should not exceed 6mm maximum.
- The amount of alloy should be adapted to the type of cylinder being used.
- The cylinder should not show any cracks.
- The size of the crucible lip should be compatible.
- The cast cone should be sufficient.

Because of the high temperature of the cylinder whilst it is being handled (on coming out of the oven) it is strongly recommended that you use a gripper and suitable gloves.

Fusing alloys is a delicate operation and the operator should devote complete attention to it. Because of this it is essential to respect some rules:

- Never add any material during the fusing process. There is a risk of hot metal being projected out.
- Observing the fusion process should not be done with the cover open and your face over the crucible.  
Monitoring is perfectly adequate using the heat-treated porthole window specially designed for this purpose.
- Use the quartz baton (and only the quartz baton) to push around the unmelted slugs and to homogenize the mixture.

**PLEASE NOTE!**

This operation requires operators to take some precautions:

- Use gloves suitable for use in heat.
- Never put your face over the crucible when doing this mixing action.
- Quickly close the cover once this operation is complete.

- The relatively high temperature of the setup at the end of the cycle can cause burns. It is essential to handle the crucible using suitable gloves and/or grippers.

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