UGICAST SPEED PLUS

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UGICAST SPEED PLUS

In the latest development from Ugin dentaire, this silica-phosphate bound carbon free investment can be used in both shock heat technique and conventional cold furnace technique with an high thermal stability. Ugicast speed plus as also the indication to be used with pressing ceramic. Due to the micrograin structure it gives clean brightly surface easy to devest and ffers a perfect fit.

MIXING RATIOS

	CONCENTRATION	RATIO LIQUID/WATER For 150g	RATIO LIQUID/WATER For 100g
NON PRECIOUS ALLOYS	100%	36ml liquid	24 ml liquid
PRECIOUS ALLOYS	50 to 60%	18 ml liquid+ 18 ml water	12 ml liquid +12ml water
PRESSABLE CERAMIC	80 %	28,8 ml liquid + 7 ,2ml water	19,2 ml liquid=4,8 ml water

All figures are given as a guide and may be varied accordind to individual technique.

When diluting Ugicast speed plus liquid use only distilled water.

EXPANSION PROPERTIES

LIQUID	SETTING EXPANSION	THERMAL EXPANSION	TOTAL EXPANSION
100%	1,2%	1,2%	2,4%
80%	0,7%	1,2%	1,9%

COLD FURNACE TECHNIQUE

Place the ring in a preheating furnace at room temperature or at maximum 150 °C.Raise to 270°C at 5/7°C per mn and hold the temperature for 20 mn.Then up to final temperature accordindg to alloy or ceramic manufacturer (max 900°C)at 5/7°C per mn and hold for 30 to 45 mn accordind to ring size.

HOT FURNACE TECHNIQUE

12 to 15 mn after powder/liquid first contact put directly ring in hot furnace at max 850°C with metal ring or at max 800°C without metal ring. Hold this temperature for 15 mn and increase at 5/7°C to the final temperature (see alloy or ceramic manufacturer instructions) (max 900°C).Holding time 45 to 60 mn depending of ring size.

DEVESTING

Cool slowly,do not dip in water !

SAFETY TYPE :

Investments contain quartz-wear adequate protection - do not breath investment dust!!!

If the **UGICAST SPEED +** Liquid comes into contact with skin wash off with soap and water. If **UGICAST SPEED +** liquid comes into contact with eyes wash out immediatly and seek medical attention !!! **UGICAST SPEED +** Liquid is alkaline.

Consult **UGICAST SPEED +** safety data sheet for further information.

All figures are approximate. When diluting **Liquid**, than use only distilled water. The above are given as a guide only and may be varied according to individual technique.

Caution : Do not bring UGIVEST SPEED + with gypsum-based materials in connection.

Storage : can be stored in a dry and frost-free room for up to 8 months.

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