# UGIVEST SPEED PLUS

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# UGIVEST SPEED PLUS

In the latest development from Ugin dentaire, this new silica phosphate bound carbon free investment developed for removable denture can be used in both shock heat technique and conventional cold furnace technique. Ugivest speed plus is suitable for all dental alloys. Due to the micro fine grain structure Ugivest speed plus gives bright castings with perfect fit and easy to deinvest.

## MIXING RATIOS

	CONCENTRATION	RATIO LIQUID and WATER For 200g	RATIO LIQUID and WATER For 100g
MODEL	55% 75%	34,5 ml liquid + 11,5 ml water 36 ml liquid + 12 ml water	17,25ml liquid + 5,75 ml water 18 ml liquid +6 ml water
MUFFLE	30% 40 %	14,4 ml liquid + 33,6 ml water 19, 2 ml liquid + 28,8 ml water	7,2 ml liquid + 16,8 ml water 9,6 ml liquid +14,4 ml water

All figures are approximate and given as a guide. The above may be varied according to individual technique.

## EXPANSION PROPERTIES

LIQUID CONCENTRATION	SETTING EXPANSION	THERMAL EXPANSION	TOTAL EXPANSION
100 %	1,6 to 1,8 %	1,2 %	2,8 to 3 %
75 %	0,7 to 1 %	1,2 %	1,9 to 2,2 %

## DUPLICATE

With silicone : let harden for 20 mn before starting wax-up. With gel : let harden 35 mn, dry model in a preheating furnace for 10 mn at 200/220°C and it in a liquid hardener (Duraform Ugin Dentaire).

## COLD FURNACE TECHNIQUE

Place the muffle into a cold furnace or at max 150°C and raise at 5/7°C to 270 °C for 20 mn, then increase at 5/7 °C to final temperature according to alloy manufacturer instructions (max 980°C) and hold at this final temperature for 30 to 45 mn depending of ring size before casting.

## HOT FURNACE TECHNIQUE

10 to 12 mn after mixing powder and liquid put the muffle in hot furnace at max 820°C without metal ring or max 850°C with metal ring, hold this temperature for 15 mn and raise at final temperature (max 980°C).

## COOLING

Allow the ring to cool slowly, do not quench in water !

## SAFETY TYPE :

Investments contain quartz-wear adequate protection - do not breath investment dust!!!

If the **UGIVEST SPEED +** Liquid comes into contact with skin wash off with soap and water. If **UGIVEST SPEED +** liquid comes into contact with eyes wash out immediately and seek medical attention !!! **UGIVEST SPEED +** Liquid is alkaline.

Consult **UGIVEST SPEED +** safety data sheet for further information.

All figures are approximate. When diluting **Liquid**, than use only distilled water. The above are given as a guide only and may be varied according to individual technique.

**Caution : Do not bring UGIVEST SPEED + with gypsum-based materials in connection. Storage : can be stored in a dry and frost-free room for up to 8 months**

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